

# VNUS® Tumescant Infiltration Pump

*Speed up your endovenous ablation procedure*



***Hands free infiltration*** – continuous flow controlled by a foot pedal

***Improve patient comfort*** – fewer needle sticks

***Increase speed*** – no need to swap needles or refill syringes

***Your ClosureFAST™ Catheter case just got FASTER***

# VNUS® Tumescence Infiltration Pump

The VNUS Tumescence Infiltration Pump facilitates rapid administration of tumescence anesthesia, thereby reducing procedure time while maximizing procedural simplicity.

The system consists of a peristaltic pump, foot pedal and single use tubing set. The flow rate can be adjusted on the Infiltration Pump to help the physician manage the infiltration speed. The foot pedal gives the physician a one handed technique controlling the fluid being infiltrated without requiring additional resources.

Sterile, single use tubings are available for use with the VNUS Tumescence Infiltration Pump. These tubing sets are individually wrapped and are sold in quantities of 10 tubing sets per box.

- **Physician comfort – eliminate repetitive action of maneuvering a syringe**
- **Variable control of flow rate on pump unit**
- **Two types of foot pedals**
  - The included Gradual Control Foot Pedal allows for hands free alteration of the flow rate during infiltration.
  - An optional Toggle Foot Pedal that provides a simple on/off mode is available



## VNUS Tumescence Infiltration Pump Ordering Information

Order#	Description
TPMP-01	VNUS Tumescence Infiltration Pump
TTUB-4	Tubing Set; sterile, 4m length(10 sets/box)
TPED-G	Additional Gradual Control Foot Pedal
TPED-T	Optional Toggle Foot Pedal

➤ **To order, call toll free: 1-888-797-VEIN (8346)**



VNUS Medical Technologies, Inc.  
5799 Fontanoso Way  
San Jose, CA 95138 USA  
Tel: 408-360-7200  
Fax: 408-365-8480  
www.vnus.com